




## **Using HiPRES to Record Your Lateral Flow Test Result**

**From 11 April 2022, you will be able to use the HiPRES app to record your lateral flow test result. Please use this guide and follow the steps to do this.**

## **Downloading the App**



You may already have the HiPRES app from taking part in saliva testing. If not you can download the latest version of the app by going to the app/play store on your device and searching for 'HiPRES.'

**If you are using Apple please ensure that you are on V11 or above of the operating system. This is the minimum operating system that the app needs to function.**

**Unfortunately the app cannot be used on Huawei devices.**

## **Registration**

If you have been using the HiPRES mobile app to take part in saliva testing, from 11 April you will be able to use the app to record the result of your lateral flow test result. You don't need to register again.

**If you haven't used the app before you need to download it (see above) then register by following the steps below:**


- Click on the blue link that says 'sign up' then enter your work email address and a password. This needs to be 8 characters long and contain one special character (do not use £), one upper case letter and one number. Then click register.
- You will receive a verification code by email. Enter this code into the box and then click verify. If you can't see this in your inbox check your junk folder. You then need to add your personal details. You can click the blue button to look up your NHS number. Or if you know what it is, type it in and click verify.
- You can find your NHS number on hospital letters, prescriptions and if you use an app for GP services. You can also look it up at [www.nhs.uk/nhs-services/online-services/find-nhs-number/](https://www.nhs.uk/nhs-services/online-services/find-nhs-number/)
- Once verified, you will see a green tick and you can click 'next' at the bottom of the screen. Then enter your phone number(s) and click next.
- You will then be asked to enter details about where you work including your manager's details. Complete all the fields and ensure contact details are accurate. Then click 'Next'.
- You will then be asked to select if you are testing yourself only (personal) or if you want to test others too (clinical and personal.) Select an option, then click next.
- Enter your email address and password. Then click 'sign in.' Select if you are testing yourself (personal testing) or someone else (clinical session.)
- If you select 'clinical and personal' you will be asked to accept some terms and conditions. Tick the box then click next.

## Recording Your Result

Once registered, login to the app and select 'LFA Pathway.' You then have the option of 'Test Self' or 'Test Others.' *(If you are a clinician, the app gives you the option to 'administer to someone else.' This means that you can submit a test result for a colleague. Please refer to the steps later in the guide if you want to use this function.)*

## Test-Self

Select the "Administer to Myself" option from the menu.



H:PRES™ Powered By IMAGEMOVER

LFA Nasal Test


Select Test Subject

Administer to Myself

Administer to Someone Else

Cancel without saving

1. Your profile is retrieved and your details displayed as shown below. Click next.



H:PRES™ Powered By IMAGEMOVER

LFA Nasal Test

Patient Details

Name: MISS Cyrilla Readett

Date of Birth: 26 Jul 2014

NHS Number: 946 434 1807

Home phone:

Mobile phone: 07770839528

Email: john.colyer@bridgeheadsoftware.com

Address:

Andover Court  
Hampton Road  
Worcester Park  
Surrey  
KT4 8EX

Ethnicity: White - British

Occupation: A2, Dept for International Development

Modify Profile...

Next >

Cancel without saving

2. You will be asked to give your consent as shown below. Click next.

**Hi-PRES™** Powered By **IMAGE MOVER**

LFA Nasal Test

**Consent**

Do you understand that COVID is a notifiable disease and that relevant authorities will be notified of a positive test?

☒ Yes

Do you consent to sharing medical data captured via this app with clinical research projects?

☒ Yes ☐ No

< Previous Next >

Cancel without saving

3. Select the kit you are using from the drop-down menu or by scanning the package of the test kit to identify the Lot Number. Alternatively, the Lot Number may be entered manually in the field below. The expiration date of the kit must also be entered. Then click next.

**Hi-PRES™** Powered By **IMAGE MOVER**

LFA Nasal Test

**Select Test Kit**

Innova SARS-CoV-2 Antigen Rapid Test ▼

Click to Scan Text

Lot number: 1344

Expiration date: 14 Jul 2022

Specimen: Nasal swab

SARS-CoV-2 Antigen Test Cartridge 1 Test

LOT X2011003

2020.11.06

2022.11.05

INNOVA CE

< Previous Next >

Cancel without saving

- Further details relating to the test then need to be collected (the cassette code from the test kit and the test location.) The cassette code is collected either manually or by scanning the barcode on the test kit. You are now ready to do the test. At this point a timer may be started to remind you when to read the result from the test. This is selected by pressing the “Start Timer” button. Alternatively, if you do not wish to use the time, you may select the “Next>” button at the appropriate time.

BridgeHead  
**H:PRES**™ Powered By **IMAGEMOVER**

LFA Nasal Test

**Process Test**

Kit: **Innova SARS-CoV-2 Antigen Rapid T**

Scan Cassette QR Code

Cassette number: Fhhbgg

Test Location: ☐ Home ☐ Work specified ▼

Please begin to administer the test by following the test instructions.

Please then wait 20 minutes before proceeding to enter the test results. You

< Previous      Next >

Cancel without saving



- If the timer is selected, a “count-down” is shown on the screen. A notification message will appear when it is time for you to read your result. When the timer expires another notification pops up. Select next to record your result.

BridgeHead  
**H:PRES**™ Powered By **IMAGEMOVER**

LFA Nasal Test

**Process Test**

Scan Cassette QR Code

Cassette number: B91672455

Test Location: Central Lancashire Area ▼

Please begin to administer the test by following the test instructions.

Please then wait 20 minutes before proceeding to enter the test results. You can set a timer to remind you when the test is done, however please note that if you choose to do this you will not be able to proceed to the next screen until the timer completes

Running... 25

< Previous      Next >

Cancel without saving

6. You will be given the option of having a copy of your result sent to you by email. Choose Yes or No and select the “Next>” button to continue.

The screenshot shows the 'Hi-PRES™ Powered By IMAGEMOVER' interface. Below the logo, it says 'LFA Nasal Test' and 'Results Communication'. A question asks 'Receive a copy of the results via email?' with two radio button options: 'Yes' (which is selected) and 'No'. At the bottom, there are three buttons: '< Previous', 'Next >', and 'Cancel without saving'.

7. You then enter the result of the test on this screen, selecting one of the three options (Positive / Negative / Invalid) then select the “Next>” button.

The screenshot shows the 'Hi-PRES™ Powered By IMAGEMOVER' interface for the 'Innova SARS-CoV-2 Antigen Rapid Test'. It displays the patient's name 'Readett, Cyrilla' and 'DOB: 26 Jul 2014'. There are three test result options, each with a visual representation of the test strip and a corresponding button: 'Positive' (with two red lines), 'Negative' (with one red line, highlighted by a green box), and 'Invalid' (with no lines). At the bottom, there are three buttons: '< Previous', 'Next >', and 'Cancel without saving'.

8. On the next screen, you are prompted to take a photograph of the test kit so that the result can be verified. When this has been done, the picture appears on the app as shown below. Click next.

The screenshot shows the HiPRES app interface. At the top, the logo 'HiPRES™ Powered By IMAGEMOVER' is displayed. Below the logo, the user's name 'Readett, Cyrilla' and date of birth 'DOB: 26 Jul 2014' are shown. A prompt 'Please take a photo of the test result.' is followed by a 'Remove' button. A photograph of a green test kit on a wooden surface is displayed. At the bottom, there are navigation buttons: '< Previous', 'Next >', and 'Cancel without saving'.

9. A summary screen is displayed. This gives brief details about you, the test kit and the result, and includes the photograph taken at the previous step. Click the “Save and Exit Session.” The test result is recorded in the HiPRES database and you are returned to the HiPRES Main Menu home page.

The screenshot shows the HiPRES app summary screen. At the top, the logo 'HiPRES™ Powered By IMAGEMOVER' is displayed. Below the logo, the word 'Summary' is shown. The user's name 'Readett, Cyrilla' is displayed, followed by the result 'Negative' in green. Below this, the test kit details are listed: 'Kit: Innova SARS-CoV-2 Antigen Rapid Test', 'Lot Number: 1344', and 'Expiration: 14 Jul 2022'. A photograph of the green test kit on a wooden surface is displayed. Below the photo, there is a 'Drag up/down to scroll' instruction. At the bottom, there are navigation buttons: '< Previous', 'Next >', and a large blue button labeled 'Save and Exit Session'. A 'Cancel without saving' link is also present at the bottom right.

## Test Others

1. Select “Administer to Others” from the menu. The following screen is displayed to collect the details of the person being tested. Enter their NHS number and date of birth and select the “Lookup” button. The HiPRES system performs a lookup on the NHS Spine.

(NHS numbers can be found at [www.nhs.uk/nhs-services/online-services/find-nhs-number/](http://www.nhs.uk/nhs-services/online-services/find-nhs-number/))

The screenshot shows the 'LFA Nasal Test' interface. At the top, it says 'HiPRES™ Powered By IMAGEMOVER'. Below this, the title 'LFA Nasal Test' is displayed. The main section is titled 'Enter NHS number or scan QR code' and contains a text input field, a 'Scan' button, and a 'Lookup' button. Below the 'Lookup' button is a 'Search by name...' input field. At the bottom right, there is a 'Cancel without saving' link.

2. If the person is successfully found, their details are displayed on the next screen. Select next.

The screenshot shows the 'LFA Nasal Test' interface with patient details displayed. At the top, it says 'HiPRES™ Powered By IMAGEMOVER'. Below this, the title 'LFA Nasal Test' is displayed. The main section is titled 'Patient Details' and contains the following information: Name: Mr Maxwell Lucey-Butler, Date of Birth: 17 Dec 2006, NHS Number: 944 930 5749, and Address: Manor Cottage, Horton Lane, Epsom, Surrey, KT19 8PF. At the bottom right, there is a 'Next >' button and a 'Cancel without saving' link.



3. You can enter additional details here about them including contact details, ethnicity and occupation. Once this information is added click next.

The screenshot shows the 'Patient Details' section of the Hi-PRES form. At the top, the Hi-PRES logo is displayed next to the text 'Powered By IMAGEMOVER'. Below this, the text 'LFA Nasal Test' is visible. The 'Patient Details' section includes the following information: Name: Mr Maxwell Lucey-Butler, Date of Birth: 17 Dec 2006, NHS Number: 944 930 5749, Email: john\_colyer@btinternet.com, Mobile: +44 7503 157350, Ethnicity: White - British, and Occupation: Accountant. Each field has a 'Change...' link next to it. At the bottom right, there is a 'Next >' button and a 'Cancel without saving' link.

Hi-PRES™ Powered By IMAGEMOVER

LFA Nasal Test

**Patient Details**

Name: Mr Maxwell Lucey-Butler

Date of Birth: 17 Dec 2006

NHS Number: 944 930 5749

Email: john\_colyer@btinternet.com

Mobile: +44 7503 157350

Ethnicity: White - British  
Change...

Occupation: Accountant  
Change...

Next >

Cancel without saving

4. Here you can enter details relating to their employment. Then the rest of the steps are the same as 'test self' (from step 2.)

The screenshot shows the 'Employment Contact Information' section of the Hi-PRES form. At the top, the Hi-PRES logo is displayed next to the text 'Powered By IMAGEMOVER'. Below this, the text 'Mr Maxwell Lucey-Butler' is visible. The 'Employment Contact Information' section includes the following information: Organisation: LSCFT Staff Member, Department: 351 L4F4 Chief Operating Officer, Sub-Department: 351 L6F7 Director Strategic Development, and Location: 351 L7F9 L9400 Strategic Development. Each field has a 'Change...' link next to it. At the bottom left, there is a '< Previous' button, and at the bottom right, there is a 'Next >' button and a 'Cancel without saving' link.

Hi-PRES™ Powered By IMAGEMOVER

Mr Maxwell Lucey-Butler

**Employment Contact Information**

Organisation: LSCFT Staff Member  
Change...

Department: 351 L4F4 Chief Operating Officer  
Change...

Sub-Department: 351 L6F7 Director Strategic Development  
Change...

Location: 351 L7F9 L9400 Strategic Development  
Change...

< Previous

Next >

Cancel without saving